

REF 6001

Jastučić za poliranje

Namjena

Jastučić za poliranje, dizajniran za uporabu na krutim endoskopima, ima površinu za poliranje koja pomaže smanjiti ostatke s leće.

Upozorenja

- Koristite samo s krutim endoskopima.
- NEMOJTE koristiti na vrlo vrućim površinama.
- NEMOJTE koristiti s kemikalijama ili otapalima.
- NEMOJTE vršiti nikakve preinake na jastučićima.
- NEMOJTE savijati endoskop tijekom poliranja.

U slučaju incidenta

Svaki ozbiljan incident koji se dogodi u vezi s ovim uređajem mora se prijaviti proizvođaču i nadležnom tijelu države članice u kojoj se nalazi korisnik i/ili pacijent.

Zbrinjavanje

Odložite u skladu sa zdravstvenom ustanovom i komunalnim smjernicama.



Domicomed.com/eifu





Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street,
Swata, BKR 4013 Malta



Casus Switzerland GmbH, Hinterbergstrasse 49,
6312 Steinhausen, Switzerland



Domico Med-Device, LLC
14241 North Fenton Rd.
Fenton, MI 48430
SAD
1.800.750.5300
www.domicomed.com



Opcionalno

